BEWARE OF THE CERMS IN THE DUST

THE SANI-GENIC POLISH MOP

-DUSTS, CLEANS, AND SHINES WITHOUT EXERTION-

it is retained in the Mon until you have gathers up the Dust from everywhere, and out or washed when necessary.

finished Dusting. It can then be shaken Moreover, it is Sanitary and Hygienic, it

We have two sizes—

Triangular Shape, Medium Size—6/6 Heart Shape, Large Size—9/6 FREE—With evey Mop we give a 40z bottle of Polish.

D.I.C.

Sole Agents

DUNEDIN





Conking's Linseed Emulsion

IS THE BEST CURE FOR

COUGHS, COLDS and CROUP

EASILY TAKEN

CHILDREN LIKE IT